

# Sampling Glassware

2012

for

## Solid Fat Content Oxidative Stability Biodiesel Analysis

### 10mm OD Pulsed NMR Tubes

For use with Bruker MiniSpec systems.

- ▶ Results are reproducible and well within acceptable experimental limits.
- ▶ Ideal for the analysis of **Solid Fat Content** in oil and fat samples.
- ▶ Caps fit securely to allow for the use of every tray position.



Catalog Number	Description	Bottom	Length	Price
NE-10PT-6-RB	Pulsed NMR Sample Tubes	Round	6"	\$98.00 / Bx/302
NE-10PT-7-FB	Pulsed NMR Sample Tubes	Flat	7"	75.00 / Bx/220
NE-10PT-7-RB	Pulsed NMR Sample Tubes auto sampler application	Round	7"	80.50 / Bx/220
NE-10PT-8-RB	Pulsed NMR Sample Tubes	Round	8"	96.00 / Bx/237
NE-310-10	Sample Tube Caps LD Polyethylene, Red	—	—	12.10 / Pk/100 65.50 / Pk/1000

### Brinkmann Rancimat Glassware

Catalog No.	Description	OD	Length	Price
SK-250504-1C	Reaction Vessel (Round Bottom)	24mm	150mm	\$232.00 / Cs/121
SK-250504-2P	Air Inlet Tube	4mm	148mm	\$69.00 / Pk/121



Reaction  
Vessel



Air Inlet Tube

### Conductivity Vessels

- ▶ 130cc capacity
- ▶ Durable (Pyrex or equivalent) laboratory glass.
- ▶ Washable and Autoclave safe.
- ▶ Outlasts plastic versions.
- ▶ Eliminates "memory effects" (which can result using porous polymers) from one sample to the next.

Catalog No.	Description	Price
SK-110604-1C	Flat Bottom Vessel, 45mm x 90mm	\$135.00 / Cs/16



New Era Enterprises, Inc.

1-800-821-4667 • Fax 1-856-697-8727

www.newera-spectro.com • cs@newera-spectro.com