



Precision NMR Sample Tubes

Sampling Accessories

Quality and value you can rely on!®


5mm Caustic Sample Cell

5mm Caustic Sample Cell

NE-321-5

Ideal for glove box applications where otherwise hazardous samples can be prepared. The cell consists of a Teflon liner with a Teflon plug and a Screw Thread sample tube configured to match Bruker 5mm probes. Agilent/Varian versions would also be available.

Your inquiries are invited.



NE-321-5	Cell, complete	\$51.99 each
NE-321-OT	Outer Tube/cap, only	34.95 each
NE-321-TL	Teflon Liner, only	18.50 each

Not sure on how you might do an NMR study on your sample. We would be happy to discuss a custom sampling solution that may suite your special sample handling. Your inquiries are welcomed and we will investigate how we might be able to help you.